

Zein

## Chemical Properties

CAS No. : 9010-66-6

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Zein is a plant-based protein polymer used to create nanofiber mats through electrospinning. Known for its excellent cell compatibility and ease of preparation, zein is suitable for use in drug delivery systems.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481