

Lactoferrin (human)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Lactoferrin (human) is an iron-binding protein that exhibits antimicrobial properties against a variety of microorganisms, including both Gram-positive and Gram-negative bacteria, viruses, protozoa, and fungi. Additionally, Lactoferrin (human) can facilitate the cultivation of Streptococcus on iron-restricted media.
Targets(IC50)	Apoptosis,Antibacterial,Antifungal

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481