

## Human Serum Albumin-AF488

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Human Serum Albumin-AF488 is a biochemical reagent conjugated with the fluorophore AF488 (Ex=488 nm; Em=496 nm). It is utilized for labeling or detecting specific antigens, molecules, or cellular structures.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481