

Carbomer 941

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Store at RT

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Carbomer 941 (CBM 941) is an acrylic polymer used as a thickening agent. It can form high-viscosity gels that provide a stable gel matrix, exhibiting excellent resistance to heat, light, and microbial contamination. Carbomer 941 aids in drug release and topical application, making it valuable in the pharmaceutical industry.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481