

Carbomer 934

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Carbomer 934 (CBM 934) is an acrylic polymer utilized as a thickening agent. It forms high-viscosity gels, providing a stable gel matrix, and demonstrates excellent resistance to heat, light, and microbial contamination. Carbomer 934 aids in drug release and topical applications, making it valuable in the pharmaceutical field.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481