

Hydroxy-PEG4-O-Boc

Chemical Properties

CAS No. : 1807530-05-7


Formula: C16H32O8

Molecular Weight: 352.42

Storage: Keep away from direct sunlight

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Hydroxy-PEG4-O-Boc is a PEG-based linker for PROTACs which joins two essential ligands, crucial for forming PROTAC molecules. This linker enables selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands linked together: one ligand binds to an E3 ubiquitin ligase, and the other targets a specific protein. They utilize the intracellular ubiquitin-proteasome system to selectively degrade these target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8375 mL	14.1876 mL	28.3752 mL
5 mM	0.5675 mL	2.8375 mL	5.675 mL
10 mM	0.2838 mL	1.4188 mL	2.8375 mL
50 mM	0.0568 mL	0.2838 mL	0.5675 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481