

UBTR008295a

Chemical Properties


CAS No. : 6947-81-5

Formula: C₁₃H₁₆O₃

Molecular Weight: 220.26

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	UBTR008295a is USP5 Inhibitor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5401 mL	22.7004 mL	45.4009 mL
5 mM	0.908 mL	4.5401 mL	9.0802 mL
10 mM	0.454 mL	2.270 mL	4.5401 mL
50 mM	0.0908 mL	0.454 mL	0.908 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481