


6beta-Hydroxystigmast-4-en-3-one

Chemical Properties

| | |
|-------------------|---|
| CAS No. : | 36450-02-9 |
| Formula: | C ₂₉ H ₄₈ O ₂ |
| Molecular Weight: | 428.69 |
| Storage: | Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA. |



Biological Description

| | |
|---------------|---|
| Description | 6beta-Hydroxystigmast-4-en-3-one is a natural product for research related to life sciences. The catalog number is TN3172 and the CAS number is 36450-02-9. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 2.3327 mL | 11.6634 mL | 23.3269 mL |
| 5 mM | 0.4665 mL | 2.3327 mL | 4.6654 mL |
| 10 mM | 0.2333 mL | 1.1663 mL | 2.3327 mL |
| 50 mM | 0.0467 mL | 0.2333 mL | 0.4665 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481