


UB 165

Chemical Properties

CAS No. : 200432-86-6
 Formula: C₁₃H₁₅ClN₂
 Molecular Weight: 234.72
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Subtype-selective nicotinic agonist
Targets(IC50)	Others,AChR

Solubility Information

Solubility	H ₂ O: <35.08 mg/mL, Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2604 mL	21.302 mL	42.604 mL
5 mM	0.8521 mL	4.2604 mL	8.5208 mL
10 mM	0.426 mL	2.1302 mL	4.2604 mL
50 mM	0.0852 mL	0.426 mL	0.8521 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481