


Anti-Phospho-53BP1 (Thr543) Polyclonal Antibody

Product Details

Reactivity: Human
Conjugation: Unconjugated
Molecular Weight: Theoretical: 220 kDa.
Purification: Affinity purification using immunogen.

Applications

Verified Activity: 1. Western blot analysis 53BP1 (Phospho-Thr543) using COLO205 whole cell lysates.



Application: WB
Recommended WB: 1:1000-3000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized phosphopeptide: human53BP1 around the phosphorylation site of Thr543
Antigen Species: human
Uniprot ID: Q12888
Synonyms: 53BP1 (p-Thr543);p-53BP1 (T543);p-53BP1 (Thr543);53BP1 (p-T543)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481