

Anti-ACTA1 Antibody (8J575)



Product Details

Ig Type:	IgG1
Reactivity:	Human,Rat,Mouse
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 103 kDa.
Clone:	8J575

Applications

Verified Activity:

1. Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using α -actinin Mouse mAb diluted at 1:200.
2. Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using α -actinin Mouse mAb diluted at 1:200.



Application:	IHC
Recommended	IHC: 1:100-200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein
Uniprot ID: P12814
Synonyms: BDPLT15;alpha-actinin-1;ACTN1;Actinin alpha 1;alpha Actinin 1;Actinin 1 Smooth Muscle;
FLJ40884;Actinin, alpha 1;F-Actin Cross-Linking Protein;FLJ54432;ACTN 1;alpha-actinin
cytoskeletal isoform;Non-Muscle alpha-Actinin-1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481