


Anti-ABCB5 Antibody (3W593)

Product Details

Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Actual: 138 kDa.
Clone:	3W593
Purification:	Affinity purification using immunogen.

Applications

Verified Activity: 1. Western blot analysis of 1) HeLa, 2) 293T, using TMAC-00041 diluted at 1:2,000.



Application:	WB
Recommended	WB: 1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Uniprot ID:	Q2M3G0
Synonyms:	ATP-binding cassette sub-family B member 5;P-glycoprotein ABCB5;ABCB5 P-gp;ABCB5

Research Background

ATP-binding cassette sub-family B member 5 also known as P-glycoprotein ABCB5 is a plasma membrane-spanning protein that in humans is encoded by the ABCB5 gene. ABCB5 is an ABC transporter and P-glycoprotein family member principally expressed in physiological skin and human malignant melanoma.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481