

Anti-Phospho-C-ABL1/2 (Tyr89/70) Polyclonal Antibody

Product Details

| | |
|-------------------|--|
| Ig Type: | IgG |
| Reactivity: | Human,Mouse,Rat |
| Conjugation: | Unconjugated |
| Molecular Weight: | Theoretical: 135 kDa. |
| Purification: | Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide. |

Applications

| | |
|--------------|------------|
| Application: | WB |
| Recommended | WB: 1:1000 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|---|
| Immunogen: | A synthesized phosphopeptide: human ABL1/2 around the phosphorylation site of Tyr89/Tyr70 |
| Antigen Species: | human |
| Uniprot ID: | P00519 & P42684 |
| Synonyms: | p-C-ABL1/2 (Y89/70);C-ABL1/2 (p-Y89/70);C-ABL1/2 (p-Tyr89/70);p-C-ABL1/2 (Tyr89/70) |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481