

Anti-Acetyl-Lysine Antibody (3R56)



Product Details

Conjugation:	Unconjugated
Clone:	3R56
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

Applications

Verified Activity:

1. Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Acetyl Lysine Monoclonal Antibody.
2. Western blot analysis of 1) Mouse Brain Tissue, 2) Hela, 3) Hela+TSA Treated using Acetyl Lysine Monoclonal Antibody.



Application:	IHC,WB
Recommended	WB: 1:1000-2000; IHC: 1:200-500

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: Purified Protein
Synonyms: Ac-K;Ac-Lysine;K-acetylated;Kac;Acetyl-K

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481