

## Anti-Acetyl-NF-kB p65 (Lys314/315) Antibody (6Z8)



### Product Details

Ig Type:	IgG1
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 65 kDa.
Clone:	6Z8

### Applications

#### Verified Activity:

1. Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Acetyl NF kB P65(K314/K315) Mouse mAb diluted at 1:200.
2. Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using Acetyl NF kB P65(K314/K315) Mouse mAb diluted at 1:200.



Application:	IHC
Recommended	IHC: 1:50-200

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: A synthesized acetylpeptide: NF kB P65 around the acetylation site of K314/K315  
Synonyms: NF-kB p65K314/315-acetylated;Ac-NF-kB p65 (Lys314/315);Ac-NF-kB p65 (K314/315);Acetyl-NF-kB p65 (K314/315);Ac-NF-kB p65K314/315;NF-kB p65K314/315ac

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481