

## Anti-Phospho-ACACA (Ser79) Polyclonal Antibody


### Product Details

Ig Type: IgG  
Reactivity: Human  
Conjugation: Unconjugated

### Applications

#### Verified Activity:

1. Western blot analysis of extracts from 293 cells untreated or treated with AICAR using Acetyl-CoA Carboxylase (Phospho-Ser79) Antibody TMAC-00066.



Application: WB

### Properties

Shipping: Shipping with blue ice.

### Antigen Details

Synonyms: ACACA (p-S79);ACACA (p-Ser79);p-ACACA (Ser79);p-ACACA (S79)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481