

## Anti-Caspase-3 Antibody (1Q94)



### Product Details

Reactivity:	Human,Mouse,Rat,Chicken
Conjugation:	Unconjugated
Molecular Weight:	Actual: 17 kDa.
Clone:	1Q94
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

### Applications

#### Verified Activity:

1. Western blot analysis of 1) HeLa, 2) 3T3, 3) Rat Brain Tissue using Active Caspase-3 Monoclonal Antibody.
2. Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using Active Caspase-3 Monoclonal Antibody.



Application:	IF,IHC,IP,WB
Recommended	WB: 1:500-1000; IHC: 1:100-200; IF: 1:50-200

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein  
Uniprot ID: P42574  
Synonyms: CPP32;SCA-1;Yama;CPP32B;apopain

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481