

Anti-Phospho-AKT1/2 (Thr308/309) Polyclonal Antibody


Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 60 kDa.
Purification:	The antibody was purified by immunogen affinity chromatography.

Applications

Verified Activity:

1. Western blot analysis of AKT1/2 (phospho-Thr308/309) expression in HeLa, PC3 whole cell lysates.



Application:	WB
Recommended	WB: 1:500-1000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic peptide: encompassing a sequence within the center region of human AKT1/2
Antigen Species:	human
Uniprot ID:	P31749 & P31751
Synonyms:	p-AKT1/2 (Thr308/309);AKT1/2 (p-Thr308/309);p-AKT1/2 (T308/309);AKT1/2 (p-T308/309)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481