

Anti-AMACR Antibody (7V845)


Product Details

Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 42 kDa.
Clone:	7V845
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

Applications

Verified Activity:

1. Western blot analysis of 1) HepG2, 2) Mouse Kidney, diluted at 1:1000.
2. IHC staining of mouse prostate adenocarcinoma tissue, diluted at 1:200.



Application:	IHC,WB
Recommended	WB: 1:1000; IHC: 1:200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Uniprot ID: Q9UHK6

Synonyms: Alpha-methylacyl-CoA racemase;2-methylacyl-CoA racemase;AMACR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481