

Anti-AMPK α 1 Antibody (5Z479)



Product Details

Ig Type:	IgG1
Reactivity:	Human,Rat,Mouse
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 63 kDa.
Clone:	5Z479

Applications

Verified Activity:

1. Western blot analysis of 1)Hela, 2) 293T, 3)3T3, 4) PC12 with AMPK α 1 Mouse mAb diluted at 1:2,000.
2. Immunohistochemical analysis of paraffin-embedded Mouse ColonTissue using AMPK α 1 Mouse mAb diluted at 1:200.



Application:	IHC,WB
Recommended	WB: 1:1000-2000; IHC: 1:50-100

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: A synthesized peptide
Uniprot ID: Q13131
Synonyms: AMPKa1;AMPK α 1;AMPK alpha1;AMPKalpha1;AMPK a1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481