


## Anti-NPPA Antibody (2T124)

## Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 16 kDa.
Clone:	2T124
Purification:	Affinity-chromatography

## Applications

Verified Activity: 1. Western blot analysis of ANP expression in PC3 cell lysate.



Application:	IP, WB
Recommended	WB: 1:1000-2000; IP: 1:20-50

## Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	A synthesized peptide: human Natriuretic peptides A
Antigen Species:	human
Uniprot ID:	P01160
Synonyms:	Prepronatriodilatin; PND; CDD ANF; Cardiodilatin; Nppa; Atrial natriuretic peptide; ANP; ANF; Atrial natriuretic factor; CDD-ANF; CDP; Natriuretic peptides A

## Research Background

Hormone that plays a key role in mediating cardio-renal homeostasis, and is involved in vascular remodeling and regulating energy metabolism.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481