

## Anti-Acetyl-HNF4A (Lys633) Polyclonal Antibody


## Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 100 kDa.
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.

## Applications

## Verified Activity:

1. Western blot analysis of mouse-lung lysate, antibody was diluted at 500. Secondary antibody was diluted at 1:20000.



Application:	WB
Recommended	WB: 1:500-2000

## Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	A synthesized acetylpeptide: human protein around the acetylation site of 633
Antigen Species:	human
Synonyms:	Ac-HNF4A (K633);Ac-HNF4A (Lys633);Acetyl-HNF4A (K633);Ac-HNF4AK633;HNF4AK633ac;HNF4AK633-acetylated

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481