

## Anti-Phospho-A-Raf (Ser299) Polyclonal Antibody

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

### Applications

Application:	IHC-P,IF(paraffin section),ELISA
Recommended	IHC-P: 1:100-300; ELISA: 1:5000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized phosphopeptide: human A-RAF around the phosphorylation site of Ser299. AA range:266-315
Antigen Species:	human
Uniprot ID:	P10398
Synonyms:	p-A-Raf (S299);p-A-Raf (Ser299);A-Raf (p-S299);A-Raf (p-Ser299)

### Research Background

A-Raf proto-oncogene, serine/threonine kinase(ARAF) Homo sapiens This proto-oncogene belongs to the RAF subfamily of the Ser/Thr protein kinase family, and maybe involved in cell growth and development. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jan 2012],

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481