

## Anti-Phospho-ATG14 (Ser29) Polyclonal Antibody

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Molecular Weight:	Actual: 55 kDa.
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.

### Applications

Application:	WB
Recommended	WB: 1:1000-2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized phosphopeptide: human Atg14 around the phosphorylation site of Ser29
Antigen Species:	Human
Uniprot ID:	Q6ZNE5
Synonyms:	p-ATG14 (S29);ATG14 (p-S29);ATG14 (p-Ser29);p-ATG14 (Ser29)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481