

## Anti-ATG7 Antibody (5A794)


### Product Details

Ig Type:	IgG1
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 78 kDa.
Clone:	5A794

### Applications

#### Verified Activity:

1. Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using ATG7 Mouse mAb diluted at 1:200.



Application:	IHC
Recommended	IHC: 1:50-200

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: ATG7
Uniprot ID:	O95352
Synonyms:	Apg 7;ATG7 autophagy related 7 homolog(S. cerevisiae);GSA7;hAGP7;APG7, S. cerevisiae, homolog of;Autophagy-related 7(yeast);Autophagy related protein 7;Ubiquitin activating enzyme E1 like protein;Atg7l;APG7L;APG7 like;APG7 autophagy 7-like(S. cerevisiae);APG7 autophagy 7 like;ATG7 autophagy related 7 homolog;GSA 7;DKFZp434N0735;Autophagy 7, S. cerevisiae, homolog of;ATG 7

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481