


Anti-Avi Tag Antibody (3Q996)

Product Details

Conjugation: Unconjugated
Clone: 3Q996
Purification: Affinity purification using immunogen.

Applications

Verified Activity: 1. Western blot analysis of Avi-Recombinant protein, using TMAC-00352 diluted at 1) 1:5,000 2) 1:10,000.



Application: WB
Recommended WB: 1:5000-10000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Synonyms: Avi Tag;GLNDIFEAQKIEWHE;AviTag;Avi

Research Background

The Avi tag is a biotin-acceptor peptide¹, GLNDIFEAQKIEWHE. The 15-residue peptide served as a substrate mimic for biotin ligase (BirA), which usually recognizes the much larger protein domain.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481