


Anti-Avi Tag Antibody (1D540)

Product Details

Ig Type: IgG1
Conjugation: Unconjugated
Clone: 1D540

Applications

Verified Activity: 1. Western blot analysis of Avi-Recombinant protein, diluted at 1) 1:5000, 2) 1:10000.



Application: WB
Recommended WB: 1:500-2000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein
Synonyms: GLNDIFEAQKIEWHE;AviTag;Avi Tag;Avi

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481