

## Anti-Phospho-AXL (Tyr698) Polyclonal Antibody

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

### Applications

Application:	ELISA,IHC,WB
Recommended	WB: 1:500-2000; IHC: 1:100-300

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized phosphopeptide: human Axl around the phosphorylation site of Y698
Antigen Species:	Human
Uniprot ID:	P30530
Synonyms:	AXL (p-Y698);p-AXL (Tyr698);AXL (p-Tyr698);p-AXL (Y698)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481