

Anti-Phospho-Beclin-1 (Ser93, 96) Polyclonal Antibody

Product Details

Reactivity: Human, Mouse, Rat
Conjugation: Unconjugated
Molecular Weight: Theoretical: 60 kDa.
Purification: Affinity purification using immunogen.

Applications

Application: WB
Recommended WB: 1:1000-3000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized phosphopeptide: human Beclin-1 around the phosphorylation site of Ser93/96
Antigen Species: Human
Uniprot ID: Q14457
Synonyms: p-Beclin-1 (S93, 96); Beclin-1 (p-S93, 96); p-Beclin-1 (Ser93, 96); Beclin-1 (p-Ser93, 96)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481