


Anti-Phospho-Brk (Tyr447) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Actual: 50 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Applications

Verified Activity: 1. Western Blot analysis of 293 cells using Phospho-Brk (Y447) Polyclonal Antibody.



Application: ELISA,IF,WB
Recommended WB: 1:500-2000; IF: 1:200-1000; ELISA: 1:5000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized phosphopeptide: human Brk around the phosphorylation site of Y447
Antigen Species: Human
Uniprot ID: Q13882
Synonyms: p-Brk (Y447);p-Brk (Tyr447);Brk (p-Y447);Brk (p-Tyr447)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481