


Anti-Serum albumin Antibody (9U460)

Product Details

Reactivity:	Bovine
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 66 kDa.
Clone:	9U460
Purification:	ProA affinity purified

Applications

Verified Activity: 1. Western blot analysis on BSA using anti-BSA Mouse mAb.



Application:	WB
Recommended	WB: 1:500-1000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	protein
Uniprot ID:	P02769

Research Background

Albumin is a soluble, monomeric protein has molecular weight of c. 66 kDa. It is primarily as a carrier protein for steroids, fatty acids, and thyroid hormones and plays a role in stabilizing extracellular fluid volume. Serum albumin protein has numerous biochemical applications including ELISAs, immunoblots, and immunohistochemistry. It is also used as a nutrient in cell and microbial culture.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481