

Anti-C1orf175 Antibody (2W831)

Product Details

Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 145 kDa.
Clone:	2W831
Purification:	ProA affinity purified

Applications

Application:	IHC,WB
Recommended	WB: 1:500-1000; IHC: 1:50-200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein
Uniprot ID:	Q68CQ1

Research Background

C1orf175 belongs to MROH family. MROH7 (Maestro Heat Like Repeat Family Member 7) is a Protein Coding gene. GO annotations related to this gene include binding. An important paralog of this gene is MROH7-TTC4.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481