

Anti-Phospho-C-ABL/ABL1 (Tyr412) Polyclonal Antibody



Product Details

Ig Type: IgG
Reactivity: Human, Mouse, Rat
Conjugation: Unconjugated
Molecular Weight: Theoretical: 122 kDa.

Applications

Verified Activity:

1. Immunofluorescence staining of methanol-fixed HeLa cells using c-Abl (phospho-Tyr412) antibody (TMAC-00514, Red).
2. Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using c-Abl (phospho-Tyr412) antibody (TMAC-00514).



Application: IF, IHC
Recommended: IHC: 1:50-100; IF: 1:100-200

Properties

Shipping: Shipping with blue ice.

Antigen Details

Uniprot ID: P00519

Synonyms: p-C-ABL/ABL1 (Tyr412);C-ABL/ABL1 (p-Tyr412);C-ABL/ABL1 (p-Y412);p-C-ABL/ABL1 (Y412)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481