


## Anti-Phospho-CACNA1C (Ser1981) Polyclonal Antibody

## Product Details

|                   |   |
|-------------------|---|
| Ig Type:          | IgG   |
| Reactivity:       | Human,Mouse,Rat   |
| Conjugation:      | Unconjugated  |
| Molecular Weight: | Theoretical: 249 kDa.   |
| Purification:     | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |

## Applications

Verified Activity: 1. Western blot analysis of lysates from 293T IFN $\alpha$ , JK T+C.



Application: WB  
Recommended WB: 1:1000

## Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

## Antigen Details

Immunogen: A synthesized phosphopeptide: human CACNA1C around the phosphorylation site of Ser1981  
Antigen Species: Human  
Uniprot ID: Q13936  
Synonyms: p-CACNA1C (S1981);CACNA1C (p-S1981);p-CACNA1C (Ser1981);CACNA1C (p-Ser1981)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481