


Anti-Phospho-CACNA1C (Ser1981) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 249 kDa.
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.

Applications

Verified Activity: 1. Western blot analysis of lysates from 293T IFNa, JK T+C.



Application: WB
Recommended WB: 1:1000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized phosphopeptide: human CACNA1C around the phosphorylation site of Ser1981
Antigen Species: human
Uniprot ID: Q13936
Synonyms: p-CACNA1C (S1981);CACNA1C (p-S1981);p-CACNA1C (Ser1981);CACNA1C (p-Ser1981)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481