


## Anti-Phospho-CASP9 (Ser196) Polyclonal Antibody

## Product Details

|                   |   |
|-------------------|---|
| Ig Type:          | IgG   |
| Reactivity:       | Human   |
| Conjugation:      | Unconjugated  |
| Molecular Weight: | Actual: 50 kDa.   |
| Purification:     | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |

## Applications

Verified Activity: 1. Western Blot analysis of 293 cells using Phospho-Caspase-9 (S196) Polyclonal Antibody.



|              |  |
|--------------|--|
| Application: | ELISA,IHC,WB                                   |
| Recommended  | WB: 1:500-2000; IHC: 1:100-300; ELISA: 1:20000 |

## Properties

|                      |   |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping:            | Shipping with blue ice.   |

## Antigen Details

|                  |   |
|------------------|---|
| Immunogen:       | A synthesized phosphopeptide: human Caspase-9 around the phosphorylation site of S196 |
| Antigen Species: | Human   |
| Uniprot ID:      | P55211  |
| Synonyms:        | p-CASP9 (Ser196);p-CASP9 (S196);CASP9 (p-Ser196);CASP9 (p-S196)                       |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481