

Anti-CTNNB1 Antibody (2B924)



Product Details

Reactivity:	Human, Mouse, Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 92 kDa.
Clone:	2B924
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

Applications

Verified Activity:

1. Western blot analysis of 1) HeLa, 2) 293T, 3) Mouse Liver Tissue, 4) Rat Liver Tissue using Catenin- β Monoclonal Antibody.
2. Immunohistochemical analysis of paraffin-embedded Human Lung carcinoma using Catenin- β Monoclonal Antibody.



Application:	IHC, WB
Recommended	WB: 1:1000-2000; IHC: 1:200-500

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein

Uniprot ID: P35222

Synonyms: PRO2286;CTNNB;armadillo;NEDSDV;EVR7;OK/SW-cl.35;beta-catenin;MRD19

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481