


Anti-CREBBP Tag Antibody (4C535)

Product Details

Conjugation: Unconjugated
Clone: 4C535
Purification: Affinity purification using immunogen.

Applications

Verified Activity: 1. Western blot analysis of CBP recombinant protein, using TMAC-00633 diluted at 1) 1:5,000 2) 1:10,000.



Application: IP, WB
Recommended WB: 1:5000-10000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Research Background

Epitope mapping within the calmodulin binding peptide (CBP) region of the TAP (tandem affinity purification) tag of human.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481