

## Anti-LAMP1 Antibody (51843)

### Product Details

Reactivity:	Human,Monkey
Conjugation:	Unconjugated
Clone:	51843
Purification:	This antibody is purified through a protein G column.

### Applications

Application:	ICC/IF,IHC,WB
Recommended	WB: 1:500-1000; IHC: 1:50-200; ICC/IF: 1:10-50

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein: human CD107a
Antigen Species:	human
Synonyms:	lysosomal-associated membrane protein 1;LAMPA;CD107a;LGP120

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481