

Anti-Phospho-gp130/IL6ST (Ser782) Polyclonal Antibody


Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 130 kDa.
Purification:	The antibody was purified by immunogen affinity chromatography.

Applications

Verified Activity:

1. Western blot analysis CD130/gp130 (Phospho-Ser782) using nocodazole treated HeLa whole cell lysates.



Application:	WB
Recommended	WB: 1:1000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	KLH conjugated synthetic phosphopeptide: corresponding to residues surrounding S782 of human CD130 protein
Antigen Species:	human
Uniprot ID:	P40189
Synonyms:	p-gp130/IL6ST (S782);gp130/IL6ST (p-S782);gp130/IL6ST (p-Ser782);p-gp130/IL6ST (Ser782)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481