

Anti-CD15 Antibody (4N938)


Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4N938

Applications

Verified Activity:

1. Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using CD 15 Mouse mAb diluted at 1:500.



Application:	IHC
Recommended	IHC: 1:50-200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide
Uniprot ID:	P22083
Synonyms:	FUC-TIV;ELFT;Fucosyltransferase 4;SSEA-1;FCT3A;LeX;SSEA1;ELAM-1 ligand fucosyltransferase; Fucosyltransferase IV(Fuc-TIV;CD15;Galactoside 3-L-fucosyltransferase;FUTIV;alpha-(1, 3)-fucosyltransferase 4;FucT-IV);4-galactosyl-N-acetylglucosaminide 3-alpha-L-fucosyltransferase;EC:2.4.1.152

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481