


Anti-FCGR2A Antibody (9Y794)

Product Details

Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Actual: 35 kDa.
Clone:	9Y794
Purification:	Affinity purification using immunogen.

Applications

Verified Activity: 1. Western blot analysis of 1) Jurkat, 2) K562 CD16 Mouse mAb(Q32) TMAC-00664.



Application:	WB
Recommended	WB: 1:1000-2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Uniprot ID:	P12318
Synonyms:	CD32A;FcγRIIA;fcγRII-a;FCGR2A;Fc gamma RIIA;FCGR2A1;FCG2;FCGR2;FCRIIA;Fc γ RIIA

Research Background

CD16 is a low affinity Fc receptor. It is a cluster of differentiation found on the surface of natural killer cells, neutrophil polymorphonuclear leukocytes, monocytes and macrophages.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481