

Anti-CD19 Antibody (1C143)



Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	1C143
Purification:	Affinity purification

Applications

Verified Activity:

1. Immunohistochemistry analysis of paraffin-embedded human normal lymph node, showing cytoplasmic localization with DAB staining using CD19 Monoclonal.
2. Western Blot analysis using CD19 Monoclonal Antibody against CD19 recombinant protein.



Application: ELISA,IF,IHC,WB

Recommended WB: 1:500-2000; IHC: 1:200-1000; ELISA: 1:10000; IF: 1:50-200

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein: human CD19
Antigen Species: human
Gene ID: 930
Uniprot ID: P15391
Synonyms: CVID3;B4;CD19 molecule

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481