

Anti-MRC1 Antibody (8Y338)

Product Details

Ig Type:	Mouse IgG1, Kappa
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 190 kDa.
Clone:	8Y338
Purification:	Protein G

Applications

Application:	ELISA,IF,WB
Recommended	WB: 1:500-2000; IF: 1:100-500; ELISA: 1:1000-5000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human protein
Antigen Species:	human
Gene ID:	4360
Uniprot ID:	P22897
Synonyms:	Human mannose receptor (hMR);Macrophage mannose receptor 1-like protein 1;MMR;MRC1; CD206;MRC1L1;CLEC13D;CLEC13DL;C-type lectin domain family 13 member D;Macrophage mannose receptor 1;C-type lectin domain family 13 member D-like

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481