

Anti-CD41 Antibody (10322)


Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	10322

Applications

Verified Activity:

1. Immunohistochemical analysis of paraffin-embedded Human Spleen using CD41 Mouse mAb diluted at 1:200.



Application:	IHC
Recommended	IHC: 1:50-200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide
Uniprot ID:	P08514
Synonyms:	GPalpha IIb;Platelet fibrinogen receptor alpha subunit;CD41b;Platelet membrane glycoprotein IIb;ITGA 2B;HPA3;CD41 antigen;CD 41;Platelet fibrinogen receptor alpha;GPalpha IIb;Integrin alpha IIb;platelet glycoprotein IIb of IIb/IIIa complex;ITA2B;Platelet specific antigen bak; Integrin, alpha 2b(platelet glycoprotein IIb of IIb/IIIa complex, antigen CD41);GP2b;Integrin alpha 2b;Integrin alpha IIb precursor;CD41a;GTA;ITGAB;GPIIb

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481