

## Anti-CD45 Antibody (1Q842)

### Product Details

Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Actual: 147 kDa.
Clone:	1Q842
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

### Applications

Verified Activity: 1. Western blot analysis of Hela, diluted at 1:2000.



Application:	WB
Recommended	WB: 1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Uniprot ID:	P08575
Synonyms:	PTPRC;LY5;T200;L-CA;B220;LCA;EC 3.1.3.48;GP180;CD45;CD45R

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481