


## Anti-Phospho-CDK1/2 (Thr14) Polyclonal Antibody

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 34 kDa.
Purification:	Affinity purification using immunogen.

### Applications

Verified Activity: 1. Western blot analysis Cdk1/2 (Phospho-Thr14) using 293 whole cell lysates.



Application:	WB
Recommended	WB: 1:1000-3000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized phosphopeptide: human Cdk1/2 around the phosphorylation site of Thr14
Antigen Species:	Human
Uniprot ID:	P06493 & P24941
Synonyms:	CDK1/2 (p-Thr14);p-CDK1/2 (T14);CDK1/2 (p-T14);p-CDK1/2 (Thr14)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481