

Anti-Collagen I Antibody (1T56)



Product Details

Ig Type:	IgG1
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 139 kDa.
Clone:	1T56

Applications

Verified Activity:

1. Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Collagen I Mouse mAb diluted at 1:200.
2. Immunohistochemical analysis of paraffin-embedded Human Skin Carcinoma Tissue using Collagen I Mouse mAb diluted at 1:200.



Application:	IHC
Recommended	IHC: 1:50-200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: A synthesized peptide: Collagen I

Uniprot ID: P02452

Synonyms: alpha 2 type I collagen;OI4;Collagen Type 1;alpha 1 type I collagen;CO1A1;COL1A1;Collagen type I alpha 2;Type I procollagen;Collagen type I alpha 1;Collagen type I;COL1A2;Collagen I alpha 1 polypeptide;Pro alpha 1(I) collagen;Osteogenesis Imperfecta Type IV;Collagen I alpha 2 polypeptide;Collagen Of Skin Tendon And Bone

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481