


Anti-COL3A1 Antibody (7P469)

Product Details

Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 138 kDa.
Clone:	7P469
Purification:	Affinity purification using immunogen.

Applications

Verified Activity: 1. Western blot analysis of Hela cells, using TMAC-00922 diluted at 1:1,000.



Application:	IHC,WB
Recommended	WB: 1:1000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Uniprot ID:	P02461
Synonyms:	CO3A1; Ehlers Danlos syndrome type IV, autosomal dominant; type III procollagen alpha 1 chain; Collagen type III alpha 1 (Ehlers Danlos syndrome type IV autosomal dominant); Collagen III alpha 1 polypeptide; COL3A1; COL 3A1; EDS4A; Collagen alpha 1(III) chain; Fetal collagen; Collagen III alpha 1 chain precursor; Collagen type III alpha; Collagen type III alpha 1; Type III collagen

Research Background

Collagen alpha-1(III) chain is a fibrillar collagen that is found in extensible connective tissues such as skin, lung, and the vascular system, frequently in association with type I collagen.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481