

Anti-CCNB1 Antibody (5K945)



Product Details

| | |
|-------------------|----------------------|
| Ig Type: | IgG1 |
| Reactivity: | Human,Rat,Mouse |
| Conjugation: | Unconjugated |
| Molecular Weight: | Theoretical: 55 kDa. |
| Clone: | 5K945 |

Applications

Verified Activity:

1. Western blot analysis of 1) Hela Cell Lysate, 2) 3T3 Cell Lysate, 3) Rat Brain Tissue Lysate using LC3A Mouse mAb diluted at 1:1000.
2. Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using Cyclin B1 Mouse mAb diluted at 1:200.



| | |
|--------------|---------------------------------|
| Application: | IHC,WB |
| Recommended | WB: 1:1000-2000; IHC: 1:100-200 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

Immunogen: A synthesized peptide
Uniprot ID: P14635
Synonyms: CCNB1;G2/mitotic-specific cyclin-B1;CCNB

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481